# Compartis® CoCr

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## D

### Intended use/indications:

Compartis CoCr are individually manufactured crowns and bridge frames for porcelain veneering by a dental technician. Compartis CoCr are custom-made. CTE (25-500 °C): 14,3 µm/m-K / (25-600 °C): 14,6 µm/m-K.

### Contraindications:

Do not use in cases of hypersensitivity to one or more of the metals contained in the alloy.

### Warnings:

- Do not inhale dust and vapours when machining.
- Ensure suitable air extraction/ventilation at the workplace and corresponding machinery

### Side effects:

Potential side effects include allergies to the metals contained in the alloy and discomforts caused by electro-chemical effects. Systemic side effects from the alloy metals are claimed in isolated cases.

### Interactions:

Avoid occlusal and approximal contact between different alloy types.

### Technical data:

Vicker's hardness	488/485 HV10
0.2% expansion limit	938/1024 MPa
Tensile strength	1337/1409 MPa
Rupture limit	3/3 %
Density	9,1 g/cm <sup>3</sup>

### Chemical composition (mass-%):

55.2% Co; 18.4% Cr; 18.4% W; 6.0% Fe; 2.0% Al.

### Soldering:

We recommend DeguDent solder N1W for repair and connection soldering of the CoCr alloy.

### Re-use:

Compartis CoCr consists of an alloy with a special composition which can only be used by laser melting and not by casting. This material should therefore not be reused!

### Processing:

Treat the inside of the crown frames with shot peening using aluminium oxide, dry off and place carefully on the stump.

In delivered state, the bridge has been shot peened out with 250 µm aluminium oxide and 3 bar.

**Recommendation:** The inner surface of the crowns has not been machined yet: small welding beads still adhere to the inner surface which can be ground down with a rosehead burr, small pear-shaped cutter or with a diamond-coated spherical burr without applying any pressure.

**Caution:** If too much pressure is applied or too much material ground away, there is a risk that the individual fit will be too loose so that the frame sways and can no longer be corrected.

Only use cross-toothed carbide cutters for the frame.

Then please scan the frames initially in circular and incisal respectively occlusal mode. For a long-life frame structure, the values measured should be no less than the following wall/connector thicknesses:

Crown frame:	front and side teeth area: 0,3 mm
Abutment:	front and side teeth area: 0,4 mm – thicknesses of 0,35 mm can be accepted in places
Bridge connector:	front teeth area: 6 mm <sup>2</sup> side teeth area: 9 mm <sup>2</sup>

Now adapt the marginal edge exactly to the preparation limits and make any necessary corrections by sanding down:

### Reducing the collar for chamfer preparation:

In the case of a highly pronounced chamfer preparation, the slightly lengthened crown may possibly result in a narrow collar beyond the preparation limits. This must be reduced first before the frame is adapted.

### Grinding the edge of the crown:

The crown edge is slightly overcontoured and can be adapted with a cutter suitable for NEM or with a polisher and ground with a thin run-out edge.

Generally comply with the following when machining CoCr work:

- Do not inhale dust and vapours when machining.
- Ensure suitable air extraction/ventilation at the workplace and corresponding machinery

As preparation for veneering, now carefully shot peen the frame with aluminium oxide and clean with the steam jet cleaner: 250 µm, 3 bar possibly less pressure so as not to damage the edges.

GB

Oxidation before ceramic firing is not necessary, as this was already carried out in the works before delivery.

#### CoCr must always be shot peened before veneering!

**Caution:** The frame has been produced as modelled in the Cercon art software. Reworking may possibly be necessary in the area of the connectors and the basal surfaces of the pontics, if Cercon art was used without smoothing function.

The frames can now be veneered straight away, heeding the firing programs specially provided for NEM.

#### Integration / fastening of Compartis restorations:

Crowns: Conventional cementing  
Multiple-unit bridges: Conventional cementing

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