

splash-guard

REF: 5562 0400

REF: 2035



Protective device used during water-cooled grinding inside milling units

Water-cooled grinding of all-ceramic materials in milling units requires reliable protection of the mechanical aspects and the magnetic plate of the respective milling units used for the grinding task.

DeguDent splash-guard was developed to afford protection of the functional aspects of the milling unit without impairing the efficacy of the grinding process.

Used for

- Conus technology using all-ceramic primary components
- Conus-supported implantology using all-ceramic primary components
- Processing of all-ceramic stock abutments in the milling unit

Technical specifications

- splash-guard corpus made of plastics
- Glass plate made of hardened glass
- Metal base made of stainless steel
- Plane-parallel fabrication and bonding of the splash-guard corpus and glass plate
- Drain groove with tube drain hole at the lower edge of the corpus
- Suitable for all popular types of milling units

Included

- splash-guard
- Tubing
- Manual

DeguDent splash-guard is available in Germany from:

DeguDent GmbH

Phone: ++49 (0) 180 23 24 555 • Fax: ++49 (0) 180 23 24 556

Internet: www.degudent.de

Sirius Dental Innovations GmbH

Phone: ++49 (0) 40 42 91 06 23 • Fax: ++49 (0) 40 42 91 06 24

E-Mail: info@sirius-dental.com

For more information about related products, trainings, publications and events please visit

www.sirius-dental.com

www.sirius-dental-training.com

DeguDent

A Dentsply International Company