

telescope-line

sirius design

telescope-line

REF 4042



Included

- 13-part diamond cutting set
 - 80, 40, 16, 8 and 4 nm cutting diamonds
 - Packaging material
 - Manual

DeguDent telescope-line is available in Germany from:

DeguDent GmbH

Phone: ++49 (0) 180 23 24 555 • Fax: ++49 (0) 180 23 24 556

Internet: www.degudent.de

Sirius Dental Innovations GmbH

Tel.: ++49 (0) 40 42 91 06 23 • Fax: ++49 (0) 40 42 91 06 24

E-Mail: info@sirius-dental.com

For more information about related products, trainings, publications and events please visit

www.sirius-dental.com

www.sirius-dental-training.com

Diamond cutting set for telescopic processing of all-ceramic dental materials

DeguDent telescope-line is a set of long-life diamond cutting instruments for easy use with all-ceramic materials.

The increasing use of all-ceramic frameworks in dental technology requires the use of cutting diamonds specifically adapted to the properties of the framework materials, facilitating efficient daily processing routines that are not tiring on the dental technician.

DeguDent telescope-line was designed for professional requirements in dental technology. This system maximizes efficiency, tool life, and safety when working on your all-ceramic framework material.

Used for

- **cercon** smartceramics®
- Telescopic finishing and adjustment of all-ceramic framework materials

Advantages and properties

- Extremely high tool life
- Effortless reduction
- Protection of framework integrity thanks to good cooling of the all-ceramic surface
- Processing efficiency
- Practical shapes
- Practical grain sizes

Technical specifications

- Natural diamonds
- Good self-cooling properties
- Long tool life thanks to a special manufacturing process
- DeguDent telescope-line diamond cutting instruments are only available with FG shafts.

Cutter designs

- Telescope cutters 80 – 4 nm
- Interlock cutters 80 – 8 nm
- Telescope cutters for thrust force distributor 80 – 8 nm

IMPORTANT

- DeguDent telescope-line diamond cutting instruments must be water-cooled during use.

DeguDent

A Dentsply International Company